MSET Demonstration Package Calibration





Mentis Sciences Education Toolkit Vision and Development History

Located in the historic Mill District of downtown Manchester, Mentis Sciences is an engineering firm which provides advanced material design and manufacturing capabilities to Department of Defense customers. Mentis specializes in the design, development and testing of advanced composite materials with a goal of providing unique flexibility, rapid development and prototyping for various composite applications.

Mentis Sciences, Inc. was founded in 1996 by John F. Dignam, following more than thirty years of service at the Army Materials Research Lab, where he served as the Director of Missile Materials. John F. Dignam spent most of his lifetime promoting national security and developing the most effective material systems to aid in countering global threats. He founded Mentis Sciences to continue promoting innovation, expertise, and emerging materials and manufacturing technologies, that will enhance U.S. security and promote economic growth.

His legacy continues under the strong and visionary leadership of John J. Dignam, who brings unique and innovative technical expertise to solving some of the nation's most daunting engineering challenges. The core values of ethics, integrity, community service, and commitment to excellence instilled by John F. Dignam live on with John J. Dignam and the Mentis team, and are apparent in every aspect of the company's structure, personality, and operations.

Mentis Sciences Internship Program recruits local high school students in good academic standing who reside in the HUBZone area of Manchester, NH. Successful youth with good attitudes and high motivation to work and learn have come through various avenues including non-traditional avenues like the Manchester Police Athletic League, The Salvation Army, and Manchester's Office of Youth Services.

Mentis makes a serious commitment of its resources to support the internship program by providing short courses in STEM related disciplines, student engineering activities and mentoring activities. In result of these courses, Mentis Sciences started to see a gap in STEM education. Biology and Life Science concepts were often the focus of science in the classroom, technology often included a smartphone app and engineering was nonexistent. Our interns and every student deserve to be introduced to STEM concepts with tools and resources that allow them to experience concepts hands-on and in a collaborative environment.

With this vision for our students, Mentis transferred skills used in their own manufacturing facility every day and descaled the concepts and tests into one integrated unit. Mentis has developed an integrated STEM toolkit that configures to complete 40 STEM tests. With limited lab space and budgets for lab testing equipment being tight, the Mentis Sciences Engineering Toolkit (MSET) departs from the high cost limited functionality of current educational testing systems.

The MSET offers a unique view into the world of material testing and physical science. Data indicates the MSET Program increased student participation in the classroom, interest in STEM careers and opportunities for females in STEM. Students develop a deep understanding in STEM, engineering and physical science concepts.



In many ways, the internship program and new shared vision has provided Mentis employees a new sense of purpose in their work. Mentis is now expanding their vision for the MSET program, beyond their own interns and are offering the MSETs STEM educational opportunities to other schools and educational partners in their community and around the United States.

Mentis believes that every student, no matter their upbringing or education status, should have the opportunity to learn, pursue their dreams and have the high-quality resources to so. This enrichment MSET program has proven to be beneficial, providing life-changing experiences for interns, students, as well as Mentis employees. We are excited to share it with you.



Calibration Introduction

In this unit students will explore the relationship between the load cells and the strain data output due to how the output strain will correlate to a force value. Students will develop the ability to map out a best fit linear curve that will correspond to a more accurate load cell the better the fit of the line is and/or the amount of data points taken along the range of force values of the load cell. The goal of this exercise is to prepare students and apply these concepts in an inquiry-based project where they will learn how variances in a calibrated system can throw off calibration and how the calibration curve can be altered to calibrate the system again.

The materials in this section have been created and organized to assist teachers in the design of lessons that use the MSET equipment and applied inquiry-based projects that are aligned with the Next Generation Science Standards, as well as the Massachusetts Science and Technology/Engineering Curriculum Framework.

Teachers

Teachers should review the Understanding by Design unit plan with particular attention to the Essential Questions students will be expected to answer by the conclusion of the unit. Throughout the lessons and experiences, teachers should assess students' progress toward their capacity answering the essential questions. Finally, teachers should use the rubric to assess students' comprehension and application of the foundational principles associate with the lesson, experiment and materials covered in the unit.

Students

It is assumed that students participating in this unit will have experience in the following areas:

- 1. The concept of a "linear" response and how that relates to a best fit curve.
- 2. Basic mechanical skills to configure the MSET based on written and visual instruction.
- 3. Basic algebra knowledge.
- 4. Inputting different scale factors for different load cells and realizing the affect they have on the output of load in the software.
- 5. Basics of how a system is calibrated and the steps that are taken to do so.



Table of Contents

5

Understanding by Design Unit Plan	6
Supporting PowerPoint Classroom Material	8
MSET Experiment Lesson Instructions	11
Poster Overview	36
Directions for Inquiry-based Project	37
Teacher Solution Key	38
Scoring Rubric for the Inquiry-based Project	39



UbD Chart – Calibration

	Desired Results							
STANDARDS/ESTABLISHED GOALS	Transfer							
Next Generation Science Standards	Students will be able to calibrate a load cell and optical encoder to assure that equipment is performing accurately.							
	Med	aning						
HS-ETS1-4. Use a computer simulation to model the impact of a proposed solution to a complex real-world problem that has numerous criteria and constraints on the interactions Forces and Interactions: HS-PS2-1. Analyze data to support the claim that	 UNDERSTANDINGS Students will understand that 1. Load cell and optical encoders are used to measure force, as well as, rotational or linear motion. 2. Load cells measure compressive, tensile, torsional, shear and other 	 What are tensile and shear forces? Why is calibration an important part of engineering and using equipment? How is shear force measured? How is tensile force measured? 						
Newton's second law of motion.	types of forces.	visition						
Massachusetts State Standards Technology and Engineering: HS-ETS1-4. Use a computer simulation to model the impact of a proposed solution to a complex real-world problem that has numerous criteria and constraints on the interactions. HS-ETS3-4(MA). Use a model to illustrate how the forces of tension, compression, torsion, and shear affect the performance of a structure. Analyze situations that involve these forces and justify the selection of materials for the given situation based on their properties.	 Students will know The importance of calibration coefficients or "sensitivities". How to apply "sensitivities" across multiple experiments. Sensors on load cells are called strain gages. Strain gages measure resistance changes that are converted to voltage and amplified. Optical encoders measure the movement of an object. Photo sensors measure the interruption of light sources as an object rotates. 	 Students will be skilled at Calibrating the MSET to assure that the load cell and optical encoders to assure the efficacy of measurements across experiments. Select appropriate calibration techniques and apply these calibration processes in multiple experiments utilizing the MSET. Analyzing data using computer software and plotting data points to examine differences. 						



Evidence

Assessment Evidence

PERFORMANCE TASK(S):

1. Calibrate a new load cell and verify accuracy of the calibrated cell. To accomplish this the students will select a load cell and perform a calibration to generate a sensitivity of the new cell. The students will select one of the MSET experiments that involve static tests and choose the appropriate calibration with differing load cells (e.g. 5kg and 10kg).

OTHER EVIDENCE:

The essential questions will be used as an entrance/exit slip to determine growth in understanding. The following computational problem will also be used as an entrance and exit slip.

See Table 3. Sensor Scale Factors MSET experiment 01 Calibration

Stage 3 – Learning Plan

Summary of Key Learning Events and Instruction

See outline of Calibration experiment summary included.

©2015 Backward Design, LLC. Used with permission.



PowerPoint Template for Instruction – Calibration



Calibration: The objective of this assignment is to calibrate a load cell and optical encoder. Calibration coefficients will be determined for both sensors, and compared to values given by the manufacturer. The calibration coefficients, called sensitivities, will be used in other experiments.



Figure 1. Side View of a Shear Load Cell







Acquisition

Acquisition

- Students will know
- The importance of calibration coefficients or "sensitivities".
- How to apply "sensitivities" across multiple
- experiments.

 Sensors on load cells are called strain gages.

 Strain gages measure resistance changes that are
- optical encoders measure the movement of an object.
- Photo sensors measure the interruption of light sources as an object rotates.
- Calibrating the MSET to assure that the load cell and optical encoders to assure the efficacy of measurements across experiments.
 Select appropriate calibration techniques and apply these calibration processes in multiple experiments utilizing the MSET.

Students will be skilled at ...

 Analyzing data using computer software and plotting data points to examine differences.

Mini-Case Study Objective 1: Scenario: After working a job during the summer, you have saved up a With your old tires that are 22 inches in diameter you travel good amount of money. You decide to get new tires on your 6omph. After putting on your new tires that are 30 inches in diameter and without recalibrating your truck, what would be your pickup truck. Your old tires were 22 inches in diameter, and the new ones are 30 inches in diameter. Will the readings on your new speed? speedometer or tachometer be affected after the installation of the new tires? If so, how would you go about re calibrating them? **Objective 2:** Formulas: You drive by a police car while on the highway and you get $\left(\frac{Diameter_{New Tires}}{Diameter_{Old Tires}}\right)mph_{Speedow}$ mphActual = pulled over. Your speedometer said that you were going 70 mph, but the police officer tells you that you were going 90 mph. You recall that you just switched out your old 22 inch mph_{Actual})Diameter_{Old Tires} diameter tires for a new set. What is the diameter of your new $Diameter_{New Tires} = \left(\frac{mp}{mph_o}\right)$

Assessment Evidence

Students will use the MSET device to conduct an experiment to calibrate a load cell and optical encoder to assure that their tires are performing accurately when given speed traveled and diameter of their tires.

Students will calibrate a new load cell and verify accuracy of the calibrated cell. To accomplish this the students will select a load cell and perform a calibration to generate a sensitivity of the new cell. The students will select one of the MSET experiments that involve static tests and choose the appropriate calibration with differing load cells (e.g. 5kg and 1okg).

MSET Experiment - Calibration

ires?



MSET Experiment Procedure– Calibration

Technical objective

The objective of this assignment is to calibrate a load cell and optical encoder. Calibration coefficients will be determined for both sensors, and compared to values given by the manufacturer. The calibration coefficients, called sensitivities, will be used in other experiments.

Background

Load cells and optical encoders are used to measure force, and rotational, or linear motion. Load cells come in a wide variety of designs, and can be used to measure compressive, tensile, torsional, shear and other forces. The MSET uses a shear style load cell for measuring tensile and compressive forces. Shear load cells are often used in digital scales. They are used in scales because they are not subject to errors associated with the location the force is applied. Figure 1 shows the cross section of a shear load cell with four sensors attached.



Figure 1: Side View of a Shear Load Cell

The sensors most commonly used on load cells are called strain gages. Strain gages are simply a small wire that changes resistance when a load is applied. The resistance change for gages is relatively small. Because of this the signal is converted to a voltage and amplified. The conversion from resistance to voltage is done by connecting the four sensors together in a configuration shown in Figure 2, called a Wheatstone bridge.





Figure 2: Wheatstone Bridge

The Wheatstone bridge requires power to work using a voltage source "Vin +" with a return current path to "V-". The output voltage generated by the load cell is shown as Vout. If the resistance of each of the our strain sensors of the load cell are equal to each other, the output voltage will be zero. Applying a load results in an imbalance of the Wheatstone bridge. This imbalance causes a change in output voltage. This output voltage is small (millivolts), but can then be amplified making it easier to measure small changes in the load.

The relationship between load and output from the load cell is either given in units of mass per output voltage expressed in Equation #1. This scale factor "*SF*" can be converted to mass per bit expressed in Equation #2. The term "bit" in the denominator is the smallest change detectable by the computer board used to acquire data. A 10-bit board will have 2^{10} or 1024 discrete increments that can be detected.

$$SF = \frac{gms.}{volt}$$
 eq. (1)

$$SF = \frac{gms.}{bit \ count}$$
 eq. (2)

Optical encoders are a common instrument used to measure the movement of something attached to it. As an example, encoders are used to measure the rotational speed of a motor. The encoder is normally attached to the motor with a coupler that transfers rotational motion from the motor to the shaft of the encoder. Figure 3 shows the external features of an optical encoder.





Figure 3: Components of an Optical Encoder

The shaft of the encoder is attached to an internal disc with slots cut into it. An optical light source or LED (light emitting diode) is used on one side of the disc while a photo sensor is attached to the other side. As the disc rotates the light source is interrupted resulting in pulses of light sensed by the photo sensor as seen in Figure 4.



Figure 4: Disc, Diode and Photo Sensor

An electronic circuit is used to condition the pulse of light as it passes a slot on the disc. Voltage pulses are then detected by a computer board configured to count pulses, represented in Figure 5. Each pulse represents a known rotational motion of the encoder. By counting pulses, the total rotation can be calculated. The MSET uses an encoder to measure rotations of the lead screw which is directly related to the linear motion of the carriage.





Figure 5: Encoder Pulses

The scale factor is calculated in the same way for the encoder as the load cell. Instead of *grams per bit count*, the encoder scale factor is in *mm per bit count*. See equation 3

$$SF = \frac{mm}{bit \ count}$$
 eq.(3)

<u>Approach</u>

This procedure covers the calibration of the 780g load cell, 5kg load cell, 20kg load cell, 50kg load cell and the motor encoder. Upon completion of the calibration, a scale factor will be generated for each sensor. These scale factors will be referenced each time the sensors are used in a future experiment.

Experiment Setup

1. Gather the following components:





- 2. In addition to the components listed above, gather a writing instrument and notebook in which to record data.
- 3. Attach the tower to the base plate as shown in the Quick Setup Guide. **Safety shield must be used for this experiment; it has been omitted from the following illustrations for clarity purposes.**
- 4. If a load cell is attached to the carriage, remove it using the provided 4mm hex wrench





Figure 6: Load Cell vs. no Load Cell

5. Attach the two thumbscrews and spacers as shown in Figure 7.



Figure 7: Placement of Thumbscrews and Spacers

Experimental Procedure

Encoder Calibration

1. Adjust the position of the carriage until the top is level with the 2cm mark on the tower scale.





Figure 8: Alignment of Carriage

2. Click "Calibrate" followed by

Start

to launch the experiment



Figure 9: Experiment Location



MSET Experiment – Calibration

18





Figure 10. Calibration Options

4. In the notebook, create a table as shown in Table 1.

Displacement (mm)	Bit Count
0	
5	
10	
15	
20	
25	

 Table 1. Encoder Calibration Data Table

- 5. Enter 0 into the "Travel" field and click
- 6. Record the output count in the cell corresponding to 0mm in the table





Figure 11. Displacement Calibration Screen

7. Next, with the speed control set to its slowest setting, displace the carriage downward 5mm using the tower scale for reference and enter "5" in the field

and click ______. Record the output count in the table for 5mm.





Figure 12: Carriage Displacement

8. Displace the carriage another 5mm and enter "10" in the field and click

. Record the output count in the table for 10mm.

- 9. Continue in this fashion until you have displaced the carriage 25mm then return the carriage to the 2cm mark
- Save Calibration
- 11. Click "yes" to save the encoder scale factor.
- 12. Next, to verify the scale factor select "Check Scale Factor" on the main menu.



General	Materials	Physics	Simple Machine	s Thermo-Fluids		
Select Calibr Check Repla Sensc	an Option ate < Scale Facto y Results or Response	γ r				
	Start					

Figure 13: Check Scale Factor Location

13. Press Encoder.	
14. If needed, press Remove Bias to zero out the reading.	
Mentis Sciences Educational Toolkit	
Update time (msec.) Travel (mm) Output count 100 0.0 0	
Remove Bias Return to Main	

Figure 14. Encoder Scale Factor Check Screen

15. Displace the carriage a known distance using the scale on the MSET tower.



- 16. Read the number displayed in "Travel" read out and compare it to the distance on the tower scale.
- 17. Ideally you would like the scale and the readout to be within a 10% error. If you find that the error is consistently above 10% it is strongly recommended that you rerun the calibration to obtain a better calibration factor. If your readings are greater than 15% then it is vital that you rerun the calibration.

$$Error = \frac{|Readout - Tower Scale|}{Tower Scale} \times 100 \qquad \text{eq. (4)}$$

Load cell Calibration-780g

1. Attach the 780g load cell to the carriage with the arrow pointed upwards as seen in Figure 15.



Figure 15: Placement of 780g Load Cell

2. Plug the load cell into port 1 at the back of the SIM.







3. Attach the hanger assembly to the load cell as seen in Figure 17.



Figure 17: Attachment of the Hanger Bolt Assembly and Wing Nut

- 4. Adjust the height of the carriage until the top is level with the 18cm mark on the tower scale.
- 5. Place the digital scale on the thumbscrew spacers and place the hook through the ring of the hanger bolt as seen in Figure 18.





Figure 18: Setup for 780g Load Cell Calibration

- 6. Continue to thread the hanger bolt into the loadcell until most of the slack between the ring and scale hook is removed but there is no load on the loadcell
- 7. Turn on the scale.
- 8. If you are using the scale for the first time, it may need to be configured to work properly with the MSET. The display should read in units of "kg". If not, press the "unit" button until "kg" is indicated. Finally, the "lock" feature must be disabled. To do this, press and hold the "Tare" button until "L-OF" appears on the screen.





Figure 19: Digital Scale Operation

9. Click "Calibrate" to launch the calibration experiment.

Mentis Sciences Educational Toolkit
General Materials Physics Simple Machines Thermo-Fluids Select an Option Calibrate Check Scale Factor Replay Results Sensor Response Sensor Response Sensor Response
Start

Figure 20: Calibration Experiment Location

10. On the next screen click	Load Cell	to begin the load cell calibration
------------------------------	-----------	------------------------------------





Figure 21. Calibration Options

11. Verify that the 780g load cell is connected to terminal 1 and select the 780g

load cell from the menu. Press



Figure 22: Sensor Selection

12. Generate a table, as shown in Table 2 in the notebook.



Weight (g)	Bit Count
0	

Table 2. Load Cell Calibration Data Table

13. Enter a load of 0g in the field labeled "Weight" and click Record the Output Count in the cell corresponding to 0g in the table.



14. Begin applying load to the load cell by inserting the 4mm allen wrench into the hole in the hanger assembly and slowly turning the assembly until the hook on the hanging scale begins to become tight as seen in Figure 23. Collect 10 data points between 0g and 700g, pick points with a relativey even spread over the range. For each data point, wait for the scale to stabilize and enter the value displayed on the scale into the MSET program and press



. Record both the measured load in grams and the output count in the table in the notebook for each point.



Figure 23: Applying a Load to the Load Cell

- 15. Once ten points have been collected, click and save the scale factor.
- 16. Next, to verify the scale factor select "Check Scale Factor" on the main menu.



Figure 24. Check Scale Factor Location





- 19. Apply load to the load cell as before. Choose three points over the 700g range. For each load, record the digital scale reading and the readout on the MSET program.
- 20. Ideally you would like the digital scale and the readout to be within a 10% error. If you find that the error is consistently above 10% it is strongly recommended that you rerun the calibration to obtain a better calibration factor. If your readings are greater than 15% then it is vital that you rerun the calibration.

$$Error = \frac{|Readout-Digital Scale|}{Digital Scale} \times 100 \qquad eq. (5)$$

5kg Load Cell Calibration

1. Follow the same steps used for calibrating the 780g load cell for the 5kg load cell.

2. This time calibrate over a range of 0-4500g.

20kg Load Cell Calibration

- 1. If your MSET is supplied with a 20kg load cell, follow the same steps used for calibrating the 780g and 5kg load cells.
- 2. Use a range of 0-19500g.

50kg Load Cell Calibration

- 1. If your MSET is supplied with a50kg load cell, follow the same steps used for calibrating the preceding load cells.
- 2. Use a range of 0-49500g

<u>Data Analysis</u>

- 1. Retrieve the data sets that were written down for each sensor.
- 2. Plot calibration data points for the load cells and encoder on the provided space below. Draw a best fit line through the data.
- 3. Determine the scale factor for each instrument by calculating the slope of the best fit line. Do this by selecting two points on the line and using equation 5. Enter the calculated scale factor into table 1 below.

30





Figure 25: Example of Slope Between Two Points

$$Slope = \frac{y_2 - y_1}{x_2 - x_1}$$
 eq. (5)

4. Scale factors are automatically calculated by the MSET program. View the scale factors calculated by the program by navigating the "C:" drive of the PC then MSET > General > Calibrate > ScaleFactors.txt

100g Load Cell SF_100g=0.086 780g Load Cell SF_750g=0.672 5kg Load Cell SF_5kg=0.210 20kg Load Cell SF_20kg=17.83 50kg Load Cell SF_50kg=43.05 Shaft Encoder SF_Disp=263.0712 Pendulum:	grams/count C grams/count C grams/count C grams/count C grams/count C millimeters	Connect to Connect to Connect to Connect to Connect to s/count	terminal terminal terminal terminal terminal	1 of 2 of 1 of 2 of 1 of	SIM SIM SIM SIM SIM									
Pendulum: Degree span negative=-90 positive=+90														
Count span count=102 @ -90 count=666 @ +90														
Cam:														
count=200 @ minumum centimeter	s=0.0													
count=800 @ maximum centimeter	s=6.0													
Encoder1-212 counts/rev														
Encoder2=180 counts/rev														
ComPort=02					1									
					1									

Figure 26: Encoder and Load Cell Scale Factors

Enter the scale factors into Table 3. Compare the scale factors calculated in question 3 with the scale factors generated by the MSET program using equation 6.

$$\% Difference = \frac{Absolute|Calculated SF - Program SF|}{\binom{Calculated SF + Program SF}{2}} \times 100 \qquad \text{eq. (6)}$$

- 5. After zeroing a force reading the load cell displays 850 bits. Calculate the applied load using the experimental scale factor obtained from the 5 Kg load cell.
- 6. Solve for the carriage travel with a starting count of 110 at 0 mm, and a final count 800.



Table 3: Sensor Scale Factors

_	Calculated Scale Factor	Program Scale Factor	% Diff
Sensor			
Encoder (mm/count)			
780g Load Cell (g/count)			
5 kg Load Cell (g/count)			
20 kg Load Cell (if applicable) (g/count)			

Encoder Calibration Plot

-										
-										
-										
					 <u> </u>	 	 		 <u> </u>	

Counts



Load Cell Calibration Plot 1



Bit Count



Load Cell Calibration Plot 2

Bit Count



Load Cell Calibration Plot 3

Bit Count



MSET - CALIBRATION



www.mset.info



Inquiry-Based Mini Projects – Calibration

After working a job during the summer, you have saved up a good amount of money. You decide to get new tires on your pickup truck. Your old tires were 22 inches in diameter, and the new ones are 30 inches in diameter. Will the readings on your speedometer or tachometer be affected after the installation of the new tires? If so, how would you go about re calibrating them?

Question 1:

With your old tires that are 22 inches in diameter you travel 60mph. After putting on your new tires that are 30 inches in diameter and without recalibrating your truck, what would be your new speed?

Formula(s):

$$mph_{Actual} = \left(\frac{Diameter_{New \ Tires}}{Diameter_{Old \ Tires}}\right) mph_{Speedometer}$$
$$Diameter_{New \ Tires} = \left(\frac{mph_{Actual}}{mph_{Speedometer}}\right) Diameter_{Old \ Tires}$$

Question 2:

You drive by a police car while on the highway and you get pulled over. Your speedometer said that you were going 70 mph, but the police officer tells you that you were going 90 mph. You recall that you just switched out your old 22 inch diameter tires for a new set. What is the diameter of your new tires?



Teacher Solution Key - Calibration

Question 1:

Given:

 $Diameter_{Old Tires} = 22in$ $Diameter_{New Tires} = 30in$ $mph_{Speedometer} = 60.mph$

Using the given information, solve for mph_{Actual}

Answer:

$$mph_{Actual} = \left(\frac{30in}{22in}\right) 60mph = 81.8mph$$

Question 2:

Given:

 $Diameter_{Old Tires} = 22in$ $mph_{Actual} = 90.mph$ $mph_{Speedometer} = 70.mph$

Solve for *Diameter*_{New Tires}

Answer:

$$Diameter_{New Tires} = \left(\frac{90. mph}{70. mph}\right) 22in = 28.3in$$



Inquiry-Based Mini Project Rubric – Calibration

	3	2	1	0	Score
Proper Use of	Used the MSET to	Struggled with using	Even at the end of	Didn't use the	
Equipment	calibrate a load cell	the MSET and	the experiment,	MSET	
	and optical encoder	getting accurate data.	struggled with the		
	to assure that	Was able to use it	used of the MSET		
	equipment is	and calibrate the	and could not		
	performing	load cells with some	accurately calibrate		
	accurately.	assistance.	the load cells.		
Accuracy of Use of	Used all terms	May have used all of	Used some of the	Didn't use any of the	
Terminology	accurately including,	the terms but one or	terms but not all of	terms in the	
	calibration	two were not used	them, or the terms	explanation of the	
	coefficient,	accurately.	were used but not	solution	
	sensitivities, load		used accurately.		
	cells, strain gauges,				
	optical encoders, and				
	slope.				
Rationale for Solution	Provided a detailed	Provided a rationale	Provided a rationale,	Didn't provide a	
	rationale for the	for their solution, but	but their explanation	rationale for their	
	choices made in their	could only briefly	was lacking the	solution	
	solution. Explanation	explain the	proper connections		
	included how to	calibration process	between the		
	properly calibrate a	and/or make	calibration process		
	load cell and/or	connections between	and the various		
	optical encoder and	the calculations they	calculations		
	how that related to	completed and a	completed.		
	calibrating the	rationale for their			
	truck's	solution to the mini-			
	speedometer/tachom	project.			



					40
	eter. A detailed explanation of why a particular method and experiment was chosen is included along with calculations for speed and diameter.				
Use of Mathematical	Accurately	May have accurately	Attempted to	Did not complete or	
Computations	completed calculations and best fit line plotting for slope then determined percent difference between MSET calculations and manual calculations to aid in calibration. Calculated stiffness of a material using the 5kg and 10kg load cells accurately and determined the sensitivity of the 10kg load cell. Calculations were used to explain the method chosen for calibration	calculated, slope, percent difference, stiffness, and sensitivity, however could not use it to explain why a specific method was chosen for calibration.	calculate slope, percent difference, stiffness, and sensitivity, but included some miscalculations.	include proper calculations for experiment.	